



VAPOR (DRY) SHIPPERS

CX100 | CXR100-11M | CXR500-0

For the safe transport of biological samples in a vapor phase.

- Advanced Concept Adsorbent that enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid
- CXR models feature replaceable absorbent material



SPECIFICATION	CX100	CXR100-11M	CXR500-0
LN2 Capacity (L)	3.7	3.7	6.4
Static Hold Time (Days)	18	16	11
Evap. Rate (Liters/Day)	0.2	0.23	0.64
DIMENSIONS			
Necktube (in/mm)	2.78 / 70.6	3.58 / 90.9	8.5 / 216
Overall Dia. (in/mm)	8.75/222.3	9.2/233.7	15.5 / 393.7
Height (in/mm)	18.5 / 472.4	19.4 / 492.7	26.9 / 683.3
WEIGHT			
Empty (lb/kg)	10 / 4.53	11.7 / 5.3	30 / 13.6
Full Weight (lb/kg)	16.5 / 7.5	18.2 / 8.3	41.5 / 18.8
CAPACITY			
No. of 1.2 ml & 2.0 ml vials (5/cane)	85	85	500
No. of 1.2 ml & 2.0 ml vials (6/cane)	102	102	500
Number of .25 cc Straws (28/cane)	784	784	1736
No. of .25 cc Straws, Bulk (2 levels) in Goblets	1640	1640	7160
Number of 1/2 straws (10/cane)	280	280	1240
No. of .5 cc Straws, Bulk (2 levels) in Goblets	800	800	3200
Number of Canisters	1	1	N / A
Wessington Part Number	8100-5100-279	8100-5100-CXR	8100-5500-CXR

* Weights are for guidance only. Actual weights may differ due to manufacturing and filling variances
 Note: Evaporation rate and static holding time are nominal. Actual rate may be affected by the nature of the contents, atmospheric conditions, container history and manufacturing tolerances.

